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English

INSTRUCTIONS FOR USE AND PRESERVATION

Dear customer,

Thank you for your trust in this Orliman product. You have acquired a quality product with a high medicinal value. Please, read the instructions carefully. If you have any doubt, get in contact with your doctor, your orthopaedist specialist or our customer service department. Orliman appreciates your choice and wishes you a speedy recovery.

REGULATION

These products satisfy the European Medical Device Directive 93/42/CEE (RD. 1591/2009). A Risk Analysis has been carried out (UNE EN ISO 14971) minimising all the existing risks. The trials have been carried out in accordance with the European regulation UNE-EN ISO 22523 of Prosthesis and Orthosis.

SUITABLE FOR

- As an orthopaedic positioning method for sprains, luxation, glenohumeral subluxation.
- Capsular displacement.
- As a conservative treatment for acute and chronic bursitis.
- Rotator cuff repair surgery.
- Shoulder endoprosthesis implantation.
- Adhesive capsulitis.
- Entrapment syndrome.
- As a method of immobilisation after implantation of a prosthetic shoulder joint and after acromioplasty.

FITTING INSTRUCTIONS

THE PRODUCT SHOULD BE ADJUSTED BY YOUR ORTHOPAEDIC TECHNICIAN OR IN ACCORDANCE WITH DOCTOR'S ORDERS.

In order to obtain a greater degree of therapeutic effectiveness in the different pathologies and to prolong the life and utility of the product, it is fundamental that the correct size adequate for each patient or user is chosen (enclosed in the packet is a diagram with all the sizes and the equivalent in centimetres). An excess compression can produce intolerance, therefore we recommend that the compression be regulated to an optimum degree.

The selected product can be one size fits all.

To fit observe the following aspects:

A- Adjustment by the technician:

The distance from the ball hand can be adjusted by the support plate to do this, remove the Velcro fastening, then slide the plate along the guide until it is in the correct position and then secure it with Velcro fasteners.

B- Patient fitting:

3 elements:

- Ambidextrous 30° or 45° pad.
- Splint with ball.
- Forearm support.

Fit the cushion to the side of the affected arm. Each side of the cushion has an inclination of 30° or 45°. Select the desired angle.

1. Secure it with the pelvic belt.
- 2-3. Attach the forearm support to the previously adjusted hand support and fasten with the hook. Position the forearm and secure it with the forearm wrist and elbow strap.
4. Fit the shoulder and side straps joining them together through the loop.
- 5-6. Tighten the shoulder and side straps.
7. Attach the carabiner to the cushion ring.
8. Finally, check the tension and tighten all of the straps of the device to ensure a perfect fit.

The sling can be fitted to the right or left arm by inverting the shoulder and torso support straps so that the shoulder strap becomes the side strap and vice versa. The fastening buckle, which can be closed from both sides, then has to be turned over.

Both support straps have been colour coded on the closure system:

Left arm: the shoulder strap is marked blue and the chest strap orange.

Right arm: the shoulder strap is marked orange and the chest strap blue.

PRECAUTIONS

The constructive material is inflammable.

Do not expose these products to situations where they could ignite. If this were the case remove them quickly and use the necessary measures to put it out. The material used is hypoallergenic, however we cannot guarantee 100% that in some determined cases it will not produce allergies, if this occurs, remove the product and see the doctor that prescribed it.

In case of minor irritations due to the use of a cotton sock, to separate the skin from contact with the material. For small inconvenience caused by sweat, we recommend using an interface to separate the skin from contact with the tissue. If rash, irritation or swelling to remove the product and consult a doctor or prosthetist. Contraindicated in open scars with swelling, redness and accumulation heat.

RECOMMENDATION-WARNINGS

These products are to be used only by those which have been medically prescribed to do so. Persons unconnected to such prescriptions should not use them. When the product is not being used, store in its original packaging. In order to dispose of the packaging and the product, strictly comply with the legal norms of your community.

The use of the orthosis is conditioned to the recommendations of the prescribing doctor, therefore it must not be used for any other means than those prescribed. For the validity of the guarantee, the issuing establishment should comply with this sheet of instructions.

MANUFACTURE-CHARACTERISTICS

The material features used in production have been tested and approved and meet all the European quality requirements. All the products are produced with high quality materials and offer an unbeatable comfort and quality of use. All products offer restraint, stability and compression for the optimum treatment of the pathologies for which they have been designed.

RECOMMENDATIONS FOR PRESERVATION AND CLEANING

Stick the velcros to each other (if the orthosis has any), wash periodically by hand with a damp cloth and a neutral soap. To dry the orthosis, use a dry towel to absorb most of the dampness and leave to dry at room temperature. Do not hang up or iron and do not expose to direct heat sources such as stoves, heaters, radiators, direct sun light etc. During use or during washing do not use alcohols, ointments or dissolvent liquids. If the orthosis is not properly dried any detergent residues could irritate the skin and deteriorate the product.

GUARANTEES

ORLIMAN, S.L.U. Guarantees all its products, as long as these have not been manipulated or altered from their original state. Does not guarantee those products in which as a result of misuse, deficiencies or breakage of any kind, its original state has been altered. If you observe any deficiency or anomaly, communicate this immediately to the establishment from which it was obtained in order for it to be changed.

ORLIMAN, S.L.U. Appreciates your choice and wishes you a speedy recovery.

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