ORLIMAN
ORTHOPEDIC PRODUCTS

Ref.: 7104-A | 6104-A

3-Tex®

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Ref.: 7104-A | 6104-A

3-Tex®

EN The manufacturing date is included in the batch number which can be found on the packaging label as LOT, in the following way: the second and third digits represent the year of manufacture and the fourth and fifth digits represent the month.

ORLIMAN
PRODUCTOS ORTOPEDICOS
ORTHOPEDIC PRODUCTS

ORLIMAN S. L. U.
C/ Ausias March, 3 - Pol. Ind. La Pobla-L´Eliana
Apdo. de correos 49 · C.P.: 46185
La Pobla de Vallbona Valencia - España (Spain)
Tel.: +34 96 272 57 04 - Fax: +34 96 275 87 00
Tel. Exportación: +34 96 274 23 33
E-mail: orto@orliman.com · Export mail: export@orliman.com
www.orliman.com

7104-A | 6104-A

Fecha de emisión/Date of issue: 2021-09

Fecha de revisión/Revision date: 2022-02 | v.02

USE AND MAINTENANCE INSTRUCTIONS **ENGLISH**

INSTRUCTIONS FOR USE, STORAGE AND WARRANTY
Dear Customer,
Thank you very much for placing your trust in an Orliman product. Please read the instructions carefully. Keep these instructions and the packaging for future reference. If you have any questions or concerns, please contact your doctor, orthopaedic specialist or our customer service department.
ORLIMAN S.L.U. guarantees all its products as long as the original configuration has not been manipulated or altered except for the intended use as described in these instructions.
If the products are used in combination with other products, replacement parts or systems, make sure they are compatible and made by Orliman®. It does not guarantee any products with altered characteristics due to improper use, defects or breakage of any kind. The statutory regulations of the country of purchase apply. Please first contact the retailer from whom you obtained the product directly in the event of a potential claim under the warranty. If any serious incidents related to the product occur, notify Orliman S.L.U. and the corresponding competent authority in your country.
Orliman would like to thank you for choosing this product and hopes you a speedy recovery.

REGULATIONS
MD This article is defined as a class I medical device. A Risk Analysis (UNE EN ISO 14971) has been carried out, minimising the existing risks. Tests have been in accordance with European Regulation UNE-EN ISO 22523 on Prostheses and Orthoses.

INDICATIONS
Arthrosis and arthritic pain and inflammation of the joints, tendinitis, minor instabilities. Congenital or acquired instabilities. As post-surgical treatment of the knee, rehabilitation processes.

FITTING INSTRUCTIONS
For best therapeutic results considering different pathologies and to extend the useful life of the product, it is essential to choose the correct size for each patient or user. Excessive compression may be intolerable; adjusting the compression to be firm yet comfortable is recommended.
If the product must be adapted, any such adjustments must be done by an orthopaedic specialist or healthcare professional legally certified to do so who must make sure the end user or person responsible for fitting the product properly understands how it works and should be used.

When fitting the product, you must adhere to the following instructions:

Ⓐ **Kneepads with lateral springs.**
Insert kneepad ensuring it meets the triangular bone (patella) window with the kneecap placed so that the patella surrounds the external outline of the kneecap. The lateral springs should be parallel, placed on the medial and lateral line of the knee.

Ⓑ **Kneepad with polycentric / joint.**
The joints should be parallel to each other, and the mechanical shaft must meet the anatomic shaft exactly.

Ⓒ **Flexo-extension kneepad.**
The joints should remain parallel to each other, the mechanical shaft should exactly meet the anatomical shaft. In order to regulate the range of mobility of flexo-extension, before adapting regulate the degrees of flexion and extension desired, which will lift the lid of the joint. Blockage takes place at an extension of 0°, the flexo-extension can be regulated at intervals of 0-15-30-60-90 degrees, by inserting the lock in the different corresponding orifices represented on the mobility scale of the joint (close the lid). If a kneepad has been chosen that incorporates velcro pins, these must be the last to be tightened.

Instructions for the orthopaedic technician:
If it is necessary to adapt the mediolateral metallic splints to the patient's anatomy:
1-Take the orthosis off the patient and alter the joint's permitted degrees by moving the pins responsible for flexion movement to the "FULL RANGE" position.
2-Move the orthosis in the direction of flexion to remove the joint (including the splints).
3-Make the modifications and place the orthosis on the patient's leg, perfecting its shape until it comes into maximum contact with the limb. Make sure the axis of the mechanical joint coincides with that of the anatomical joint.
NOTE: Exercise care with the anatomy of the femoral condyles: ensure the upper and lower ends of the splints do not squeeze the thigh or calf. Carry out this operation on both the internal and external joints, keeping their anatomical differences in mind.

⚠ **PRECAUTIONS**

Before each use, check that all product components are present as per the fitting process. Periodically check the conditions of the product. If you observe any defect or anomaly, immediately report it to the issuing establishment.

This product is made of inflammable material. Do not expose the products to situations that could set them on fire. In the event of a fire, quickly get them off your body and use the proper resources to extinguish the fire.

To avoid minor discomfort caused by sweating, we recommend using some type of cotton fabric to separate the skin from contact with the product material. For discomfort such as chafing, irritation and swelling, remove the product and see a doctor or orthopaedic specialist. The product should only be used on healthy skin. It is not recommended for use over open scars with swelling, redness or hotspots.

Products marked with the LATEX symbol contain natural rubber latex and can cause allergic reactions in people sensitive to latex.

Products marked with the MR symbol contain ferromagnetic components and, therefore, extreme precaution must be taken if you undergo an MRI scan or are exposed to radiation associated with diagnostic or therapeutic procedures.

RECOMMENDATIONS-WARNINGS
The use of these products is conditioned by the indications. Although the product is not defined as a single-use device, using it on a single patient only is recommended and only for the intended purposes as described in these instructions or by a healthcare professional.
When disposing of the product and its packaging, you must strictly adhere to the legal regulations in your community.

RECOMMENDATIONS FOR STORAGE AND WASHING
When not using the product, store it in the original packaging in a dry place at room temperature. Stick the Velcro to each other (if the orthotic device has them), frequently wash by hand with warm water (30º C max.) and mild soap. To dry the product, use a dry towel to absorb as much moisture as possible and let it dry at room temperature. Do not hang it up or iron the product and do not expose it to direct heat sources such as stoves, dryers, direct sun exposure, etc. When using or cleaning the product, do not use abrasive or corrosive substances, alcohol, ointments or liquid solvents. If not dried off properly, the detergent residue may irritate the skin and cause the product to deteriorate.

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